

Tested Devices – Picterus Jaundice Pro

The Picterus Jaundice Pro app has been designed and built to work on any modern smart phone, the minimum device requirements for both Android and Apple smartphone are:

Android Smartphone

Minimum Requirements

- OS: ≥ Android 8.0 (Oreo)
- CPU: Octa-core or better
- RAM: ≥ 4 GB
- Internal Storage Capacity: ≥ 32 GB (UFS 2.0)
- Display: ≥ 5.1" (577 ppi)
- Main Camera (rear): ≥12 MP (jpeg/raw) and LED flash
- Connectivity: Internet Access - Wi-Fi or mobile data (≥ 2G)

Apple iPhone

Minimum Requirements

- OS: ≥ iOS 14.3
- CPU: Hexa-core or better
- RAM: ≥ 2 GB
- Internal Storage Capacity: ≥ 64 GB
- Display: ≥ 4.7" Retina (326 ppi)
- Main Camera (rear): ≥ 12 MP and LED flash
- Connectivity: Internet Access - Wi-Fi or mobile data (≥ 2G)

Despite this, there are sometimes issues with using the app on certain specific phone models.

Any phone model that is not on either of the following lists has not been tested by us and there is a chance that there might be issues using the app. If you attempt to use such a device and you experience problems we would greatly appreciate if you would let us know by contacting us at support@picterus.com.

The following devices have been tested by Picterus and have been ensured to work with the Picterus Jaundice Pro app.